COOK
MEDICAL

CE 0123

EN
3

Oasis® Extracellular Matrix

Instructions for Use

CS
5

Extracelulární matrice Oasis®

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DA
7

Oasis® Extracellulær Matrix

Brugsanvisning

DE
8

Oasis® Extrazellulären Matrix

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EL
11

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OASIS® EXTRACELLULAR MATRIX

INTENDED USE

OASIS® EXTRACELLULAR MATRIX (ECM) is indicated for the management of wounds including:

- Partial and full-thickness wounds
- Pressure ulcers
- Venous ulcers
- Chronic vascular ulcers
- Diabetic ulcers
- Trauma wounds (abrasions, lacerations, skin tears)
- Second-degree burns
- Donor sites/grafts
- Draining wounds
- Surgical wounds (post-Mohs' surgery, post-laser surgery, podiatric, wound dehiscence)

OASIS ECM is supplied sterile in peel-open packages and is intended for one-time use.

CAUTION: This product is intended for use by trained medical professionals.

COMPOSITION

OASIS ECM is derived using a process that retains the natural composition of matrix molecules such as collagen (Types I, III, IV, VI), glycosaminoglycans (hyaluronic acid, chondroitin sulfate A and B, heparin, and heparan sulfate), glycoproteins (fibronectin) and growth factors (FGF-2, TGF- β).^{1, 2, 3}

CONTRAINDICATIONS

This device is derived from a porcine source and should not be used in patients with known sensitivity to porcine material. This device is not indicated for use in third degree burns.

PRECAUTIONS

- Do not re-sterilize. Discard all open and unused portions of the OASIS ECM.
- Device is sterile if the package is dry, unopened and undamaged. Do not use if the package seal is broken.
- The device must be used prior to the expiration date.
- Discard device if mishandling has caused possible damage or contamination.
- OASIS ECM should not be applied until excessive exudate, bleeding, acute swelling, and infection is controlled.

POTENTIAL COMPLICATIONS

The following complications are possible. If any of these conditions occur, the device should be removed.

- Infection
- Chronic inflammation (Initial application of wound dressings may be associated with transient, mild, localized inflammation.)
- Allergic reaction
- Excessive redness, pain, swelling, or blistering

STORAGE

This device should be stored in a clean, dry location at room temperature.

STERILIZATION

This device has been sterilized with ethylene oxide.

SUGGESTED INSTRUCTIONS FOR USING OASIS EXTRACELLULAR MATRIX

These recommendations are designed to serve only as a general guideline. They are not intended to supersede institutional protocols or professional clinical judgment concerning patient care.

NOTE: Always handle the OASIS ECM using aseptic technique.

1. Prepare wound area using standard of care methods to ensure wound is free of debris and necrotic tissue. If necessary, surgically debride and treat the wound bed with standard of care techniques to ensure the wound shows granulation tissue and the edges contain viable tissue.
2. Cut the dry sheet to a size slightly larger than the outline of the wound area. If the wound is larger than a single sheet, then multiple sheets may be used. Overlap adjoining sheets to provide coverage of the entire wound.
3. Apply OASIS ECM in its dry state over the wound. Rehydrate the device with sterile saline or other isotonic solution. Alternatively, rehydrate OASIS ECM by placing it in a bowl of sterile saline or other isotonic solution for at least one (1) minute prior to use.
4. Place the edge of the device in contact with the intact tissue. Smooth the OASIS ECM into place to ensure the sheet is in contact with the underlying wound bed.
5. After application, use an appropriate, non-adherent, secondary dressing to maintain a moist wound environment. Change the secondary dressing as needed to maintain a moist, clean wound area. Frequency of secondary dressing change will be as needed.
6. Discard any unused portions according to institutional guidelines for biological waste. Reapply, as needed, if the OASIS ECM is no longer visible. Typically, reapplication is needed every 3-7 days until the wound is re-epithelialized.

IMPORTANT: As healing occurs, sections of the device may darken in color or form into a caramel-colored gel which may begin to peel from the area surrounding the wound. Any product that has lifted from the skin may be trimmed, taking care not to forcibly remove sections of the OASIS ECM that may adhere to the wound.

REFERENCES

1. Hodde J, Janis A, Ernst D, Zopf D, Sherman D, Johnson C. Effects of sterilization on an extracellular matrix scaffold: Part I. Composition and matrix architecture. *J Mater Sci Mater Med* 2007;18:537-543.
2. Hodde JP, Badylak SF, Brightman AO, et al. "Glycosaminoglycan Content of Small Intestinal Submucosa: A Bioscaffold for Tissue Replacement." *Tissue Engineering*, 1996; 2: 209-217.
3. Data on File, Cook Biotech, Inc.

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REFERENCES: See English "REFERENCES"

Keep dry

Manufacturer

Use-by date

Attention, see instructions for use

Do not re-use

Temperature limit

Magnetic resonance safe

STERILE

EO

Sterilized using ethylene oxide

MANUFACTURER

COOK BIOTECH
INCORPORATED
1425 Innovation Place
West Lafayette, IN 47906 U.S.A.
Phone: 812 339-2235
Toll Free: 800 457-4500
Toll Free Fax: 800 554-8335

EC REPRESENTATIVE

COOK IRELAND
O'Halloran Road
National Technological Park
Limerick, IRELAND
Phone: +353 61 334440

AUSTRALIAN SPONSOR

WILLIAM A. COOK
AUSTRALIA PTY. LTD.
Brisbane Technology Park
95 Brandl Street
Eight Mile Plains
Brisbane, QLD 4113
Australia
Phone: +61 7 38 41 11 88

www.cookmedical.com
© COOK BIOTECH 2018

September 2018